
Premium Grade Sulfo-SMCC (Pierce™)

产品编号: PG82085

保存条件: Store at -20°C.

产品介绍: Thermo Scientific Pierce Premium Grade Sulfo-SMCC is our highest quality formulation of amine-to-sulphydryl crosslinker, specially characterized for applications where product integrity and risk minimization are paramount.

Compared to the standard grade product, Premium Grade Sulfo-SMCC provides more clearly defined quality and product support by including (a) increased analytical testing and product characterization, (b) greater batch-specific information and quality assurance review, (c) extensive lot sample retention and (d) change control notification. Sulfo-SMCC is a water-soluble, non-cleavable and membrane-impermeable amine-to-sulphydryl crosslinker that contains NHS-ester and maleimide reactive groups at opposite ends of a medium-length cyclohexane spacer arm (8.3 angstroms). NHS esters react with primary amines at pH 7-9 to form stable amide bonds. Maleimides react with sulphydryl groups at pH 6.5-7.5 to form stable thioether bonds. The maleimide groups of Sulfo-SMCC are unusually stable up to pH 7.5 because of the cyclohexane bridge in the spacer arm. Because it contains the hydrophilic sulfonyl moiety, Sulfo-SMCC is soluble up to nearly 10 mM in water and many commonly used buffers, thus avoiding the use of organic solvents which may perturb protein structure.

Features of Sulfo-SMCC:

- Sulfo-SMCC—a water-soluble amine-to-sulphydryl crosslinker with a stabilizing cyclohexane spacer arm
- High quality—identity and purity confirmed by several tests, including quantitative NMR
- Product integrity—enhanced level of testing and characterization compared to standard grade
- Lot retention—ample supply of past lots retained to ensure future process testing
- Change management—Change Control Notification (CCN) service
- Consistent manufacture—batch-specific manufacturing documentation review


Applications:

- Enzyme labelling of antibodies - both enzyme and antibody specificity can be preserved
- Create specific bioconjugates via one- or two-step crosslinking reactions
- Create sulphydryl-reactive, maleimide-activated carrier proteins for coupling haptens

Product References:

Crosslinker Application Guide -- search for recent literature references for this product

Chemical structure of Sulfo-SMCC crosslinking reagent



Sulfo-SMCC

Sulfosuccinimidyl 4-(*N*-maleimidomethyl)cyclohexane-1-carboxylate

MW 436.37

Spacer Arm 8.3 Å