

bs-16471R

[Primary Antibody]

Hepatitis C Virus 2a E1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Hepatitis C Virus 2a E1

Immunogen: KLH conjugated synthetic peptide derived from Hepatitis C Virus 2a E1: 151-250/3033.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: Hepatitis C virus subtype 2a)

Predicted MW.: 21 kDa