

bs-18336R

[Primary Antibody]

LONP2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 83752</p> <p>Target: LONP2</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human LONP2: 531-630/852.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: ATP-dependent serine protease that mediates the selective degradation of misfolded and unassembled polypeptides in the peroxisomal matrix. Necessary for type 2 peroxisome targeting signal (PTS2)-containing protein processing and facilitates peroxisome matrix protein import.</p>	<p>Isotype: IgG</p> <p>SWISS: Q86WA8</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Dog)</p> <p>Predicted MW.: 22 kDa</p> <p>Subcellular Location: Cytoplasm</p>
--	--	---