

bs-19461R

[Primary Antibody]

PSMA8 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 143471</p> <p>Target: PSMA8</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PSMA8: 121-220/256.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8TAA3</p>	<p>Applications: WB (1:100-1000) IHC-P (1:200-500) IHC-F (1:200-500) IF (1:200-500) ICC/IF (1:100-500)</p> <p>Reactivity: Human, Mouse, Rat (predicted: Rabbit, Pig, Sheep, Cow, Dog, Horse, Cynomolgus Monkey, Rhesus monkey, Opossum)</p> <p>Predicted MW.: 29 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
--	--	---