

bs-14221R

[Primary Antibody]

DDHD1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 80821</p> <p>Target: DDHD1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human DDHD1: 751-850/900.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Phosphatidic acid is released following cell activation and functions as a second messenger in several signaling pathways. DDHD1 is a lipase that catalyzes degradation of phosphatidic acid and attenuates cell activation.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8NEL9</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog)</p> <p>Predicted MW.: 100 kDa</p> <p>Subcellular Location: Cytoplasm</p>
--	--	---