bs-14221R

[Primary Antibody]

DDHD1 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 80821 **SWISS: Q8NEL9**

Target: DDHD1

Immunogen: KLH conjugated synthetic peptide derived from human DDHD1:

751-850/900.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: Phosphatidic acid is released following cell activation and

functions as a second messenger in several signaling pathways. DDHD1 is a lipase that catalyzes degradation of phosphatidic acid

and attenuates cell activation.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) **IF** (1:100-500) ICC/IF (1:100-500) **ELISA** (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Rabbit, Pig, Sheep,

Cow, Dog)

Predicted MW.: 100 kDa

Subcellular Cytoplasm