

**bs-11303R****[ Primary Antibody ]****Orexin B Rabbit pAb**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**— DATASHEET —**

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                      |                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host:</b> Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | <b>Isotype:</b> IgG  | <b>Applications:</b> <b>WB</b> (1:500-2000)<br><b>IHC-P</b> (1:100-500)<br><b>IHC-F</b> (1:100-500)<br><b>IF</b> (1:100-500)<br><b>ICC/IF</b> (1:100-500)<br><b>ELISA</b> (1:5000-10000)<br><br><b>Reactivity:</b> (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog)<br><br><b>Predicted MW.:</b> 3 kDa<br><br><b>Subcellular Location:</b> Cytoplasm |
| <b>Clonality:</b> Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                      |                                                                                                                                                                                                                                                                                                                                                                  |
| <b>GeneID:</b> 3060                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | <b>SWISS:</b> O43612 |                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Target:</b> Orexin B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                      |                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Immunogen:</b> KLH conjugated synthetic peptide derived from human Orexin B: 51-130/130.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |                      |                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Purification:</b> affinity purified by Protein A                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                      |                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Concentration:</b> 1mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                      |                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                                                                                                                                                                                                                                                                                                                                                                                                                         |                      |                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Background:</b> Prepro-orexin is 130 amino acid long peptide with a putative 33 AA secretory sequence, a hydrophobic core followed by residues with small polar side chains. The expression was detected in brain and to a small extent in testis. These neuropeptides bind and activate two closely related Orexin receptors—G-protein coupled receptors (GPCRs) OX1R and OX2R. Orexins (Orexin A and Orexin B) are a family of hypothalamic neuropeptides selectively expressed in the hypothalamus. Orexin A and Orexin B are derived from the same precursor (Prepro-orexin) by proteolytic cleavage. |                      |                                                                                                                                                                                                                                                                                                                                                                  |