

bs-18369R

[Primary Antibody]

LRRC18 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 474354</p> <p>Target: LRRC18</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human LRRC18: 31-130/261.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: LRRC18 is a 261 amino acid protein that contains seven LRR repeats. Localized to the cytoplasm, LRRC18 may be involved in spermatogenesis and sperm maturation. LRRC18 is expressed as two isoforms produced by alternative splicing and is encoded by a gene mapping to human chromosome 10.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8N456</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Horse)</p> <p>Predicted MW.: 30 kDa</p> <p>Subcellular Location: Cytoplasm</p>
---	--	---