

bs-21030R

[Primary Antibody]

SLCO5A1/SLC21A15 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 81796</p> <p>Target: SLCO5A1/SLC21A15</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SLCO5A1/SLC21A15: 581-680/848. < Extracellular ></p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene encodes a 12 transmembrane domain protein that is a member of the solute carrier organic anion transporter superfamily. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2009]</p>	<p>Isotype: IgG</p> <p>SWISS: Q9H2Y9</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 92 kDa</p> <p>Subcellular Location: Cell membrane</p>
--	--	--