
SLITL2 Rabbit pAb

Catalog Number: bs-21034R

Target Protein: SLITL2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat)

Predicted MW: 72 kDa

Entrez Gene: 114990

Swiss Prot: Q6EMK4

Source: KLH conjugated synthetic peptide derived from human SLITL2: 451-550/673.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: VASN (Vasorin) is a Protein Coding gene. GO annotations related to this gene include transforming growth factor beta binding.