

---

## LANPL Rabbit pAb

Catalog Number: bs-6862R

Target Protein: LANPL

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)

Predicted MW: 30 kDa

Entrez Gene: 81611

Swiss Prot: Q9BTT0

Source: KLH conjugated synthetic peptide derived from human LANPL: 55-150/268.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Inhibits activity of protein phosphatase 2A (PP2A). Does not inhibit protein phosphatase 1.

May play a role in cerebellar development and synaptogenesis process by modulating PP2A activity.

Tissue specificity Expressed in peripheral blood leukocytes, colon, small intestine, prostate, thymus, spleen, skeletal muscle, liver and kidney.