
Thioester-containing protein 1 Rabbit pAb

Catalog Number: bs-0632R

Target Protein: Thioester-containing protein 1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Anopheles stephensi)

Predicted MW: 29 kDa

Source: KLH conjugated synthetic peptide derived from Anopheles stephensi Thioester-containing protein 1: 101-180/260.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: No Data Available.