

**bs-0632R**

**[ Primary Antibody ]**

## Thioester-containing protein 1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>Target:</b> Thioester-containing protein 1</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from Anopheles stephensi Thioester-containing protein 1: 101-180/260.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> No Data Available.</p>	<p><b>Isotype:</b> IgG</p> <p><b>Applications:</b> <b>WB</b> (1:500-2000) <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:100-500) <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Anopheles stephensi)</p> <p><b>Predicted MW.:</b> 29 kDa</p>
--	--