

bs-11188R

[Primary Antibody]

CEBP gamma Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 1054</p> <p>Target: CEBP gamma</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human CEBP gamma: 51-150/150.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: CEBP gamma is a transcription factor that binds to the enhancer element PRE-I (positive regulatory element-I) of the IL-4 gene. It might change the DNA-binding specificity of other transcription factors and recruit them to unusual DNA sites.</p>	<p>Isotype: IgG</p> <p>SWISS: P53567</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 16 kDa</p> <p>Subcellular Location: Nucleus</p>
---	--	--