

bs-12783R

[Primary Antibody]

VSTM2B Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 342865 Target: VSTM2B Immunogen: KLH conjugated synthetic peptide derived from human VSTM2B: 31-130/285. < Extracellular > Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: Contains 1 Ig-like V-type (immunoglobulin-like) domain.	Isotype: IgG SWISS: A6NLU5	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog) Predicted MW.: 27 kDa Subcellular Location: Cell membrane
---	--	--