

bs-11392R

[Primary Antibody]

AAV5 capsid VP2 protein Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: AAV5 capsid VP2 protein

Immunogen: KLH conjugated synthetic peptide derived from Adeno-Associated Virus 5 capsid protein VP1: 511-560/724.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Adenoviruses are DNA viruses that cause upper respiratory tract infections, conjunctivitis, and other infections in humans. Adenoviruses were discovered in 1953. About 47 different types have been identified since then, and about half of them are believed to cause human diseases. Infants and children are most commonly affected by adenoviruses. Adenovirus infections can occur throughout the year, but seem to be most common from fall to spring. Adenoviruses are responsible for 3-5% of acute respiratory infections in children and 2% of respiratory illnesses in adults.

Applications: **IHC-P** (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: AAV5)

Predicted MW.: 80 kDa