
DEGS2 Rabbit pAb

Catalog Number: bs-14261R

Target Protein: DEGS2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Cow, Dog, Horse)

Predicted MW: 37 kDa

Entrez Gene: 123099

Swiss Prot: Q6QHC5

Source: KLH conjugated synthetic peptide derived from human DEGS2: 2-100/323.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a bifunctional enzyme that is involved in the biosynthesis of phytosphingolipids in human skin and in other phyt