

bs-18412R

[Primary Antibody]

LRRC8D Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 55144 Target: LRRC8D Immunogen: KLH conjugated synthetic peptide derived from human LRRC8D: 44-153/858. < Extracellular > Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: LRRC8D contains 12 LRR (leucine-rich) repeats. The function of this protein is unknown.	Isotype: IgG SWISS: Q7L1W4	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit) Predicted MW.: 98 kDa Subcellular Location: Cell membrane
--	--	--