

**bs-24097R**

**[ Primary Antibody ]**

## Neuromedin U Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>Target:</b> Neuromedin U</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from mouse Neuromedin U: 144-166/174.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> Stimulates muscle contractions of specific regions of the gastrointestinal tract. In humans, NmU stimulates contractions of the ileum and urinary bladder.</p>	<p><b>Isotype:</b> IgG</p> <p><b>Applications:</b> <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:100-500) <b>ICC/IF</b> (1:100-500)</p> <p><b>Reactivity:</b> (predicted: Mouse)</p> <p><b>Predicted MW.:</b> 20 kDa</p> <p><b>Subcellular Location:</b> Secreted</p>
--	--