
Rabbit Anti-Chicken IgG H&L, Gold conjugated

Catalog Number: bs-0310R-Gold

Target Protein: Rabbit Anti-Chicken IgG H&L

Concentration: 0.4 mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Reactivity: Chicken

Purification: affinity purified by Protein A

Storage: 20mM TBS (pH8.0) with 1% BSA and 0.02% Proclin300.

Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.

Background: Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10¹¹ variants.

PRODUCT SPECIFIC PUBLICATIONS

[IF=3.939] Xiangguang Li . et al. Production and characteristics of a novel chicken egg yolk antibody (IgY) against periodontitis-associated pathogens. J Oral Microbiol. 2020;12(1):1831374 IF ; actinomycetemcomitans . 33144924