
Neutrophil/Ly-6B Rabbit pAb

Catalog Number: bs-19701R

Target Protein: Neutrophil/Ly-6B

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),
ELISA (1:5000-10000)

Reactivity: (predicted:Rat)

Predicted MW: 14 kDa

Swiss Prot: Q63317

Source: KLH conjugated synthetic peptide derived from Rat Ly-6B: 27-100/135.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.