

bs-19701R

[Primary Antibody]

Neutrophil/Ly-6B Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

SWISS: Q63317

Target: Neutrophil/Ly-6B

Immunogen: KLH conjugated synthetic peptide derived from Rat Ly-6B: 27-100/135.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: **WB** (1:500-2000)

IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Rat)

**Predicted
MW.:** 14 kDa

**Subcellular
Location:** Cell membrane