

bs-21050R

[Primary Antibody]

PP11 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|--|---|
| <p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 8909</p> <p>Target: PP11</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PP11: 161-260/410.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene encodes a protein with protease activity and is expressed in the placenta. The protein may be useful as a tumor marker. Multiple alternatively spliced transcript variants have been found for this protein. [provided by RefSeq, Feb 2010]</p> | <p>Isotype: IgG</p> <p>SWISS: P21128</p> | <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat)</p> <p>Predicted MW.: 47 kDa</p> <p>Subcellular Location: Secreted</p> |
|--|--|---|