

**bs-21050R**

**[ Primary Antibody ]**

## PP11 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|  |  |   |
|--|--|---|
| <p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 8909</p> <p><b>Target:</b> PP11</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human PP11: 161-260/410.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br/>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> This gene encodes a protein with protease activity and is expressed in the placenta. The protein may be useful as a tumor marker. Multiple alternatively spliced transcript variants have been found for this protein. [provided by RefSeq, Feb 2010]</p> | <p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> P21128</p> | <p><b>Applications:</b> <b>IHC-P</b> (1:100-500)<br/><b>IHC-F</b> (1:100-500)<br/><b>IF</b> (1:100-500)<br/><b>ICC/IF</b> (1:100-500)<br/><b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Human, Mouse, Rat)</p> <p><b>Predicted MW.:</b> 47 kDa</p> <p><b>Subcellular Location:</b> Secreted</p> |
|--|--|---|