

**bs-10298R**

**[ Primary Antibody ]**

## H7N9 NA Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** H7N9 NA

**Immunogen:** KLH conjugated synthetic peptide derived from H7N9 neuraminidase: 1-100/465.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** WB (1:500-2000)

**IHC-P** (1:100-500)

**IHC-F** (1:100-500)

**IF** (1:100-500)

**ICC/IF** (1:100-500)

**ELISA** (1:5000-10000)

**Reactivity:** (predicted: Chicken, H7N9)

**Predicted MW.:** 51 kDa

**Subcellular Location:** Cytoplasm ,Nucleus