

bs-6979R

[Primary Antibody]

hnRNP R Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 10236</p> <p>Target: hnRNP R</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human hnRNP R: 535-633/633.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Component of ribonucleosomes, which are complexes of at least 20 other different heterogenous nuclear ribonucleoproteins (hnRNP). hnRNP play an important role in processing of precursor mRNA in the nucleus.</p>	<p>Isotype: IgG</p> <p>SWISS: O43390</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:50-200) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Dog, Horse)</p> <p>Predicted MW.: 70 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
--	--	---