

**bs-0339G**

**[ Primary Antibody ]**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

## **T3 Goat pAb**

### **— DATASHEET —**

**Host:** Goat

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** T3

**Immunogen:** BSA conjugated T3: T3.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** T3 (Triiodothyroxine) although only a minority of the total thyroxine, but the T3 is very important to check thyroid function. In the biochemical function of T4 than strong, and in the case of hyperthyroidism, T3 concentration in the blood are often marked increase (or even up to 5 times the normal).

**Applications:** ELISA (1:5000-10000)

**Reactivity:** (predicted: T3)

**Predicted MW.:** 0.65097 kDa