



www.bioss.com.cn
sales@bioss.com.cn
techsupport@bioss.com.cn
400-901-9800

Human Tau Mouse mAb

Cat. No:	V5213
Name:	Human Tau Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	7A2
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant human Tau-4 protein (bs-4133P)
Pairing Information:	

Pairing Information:

ELISA	Detection																	
	V5203	V5206	V5214	V5215	V5216	V5218	V5219	V5232	V5233	V5234	V5236	V5237	W5201	V5209	V5210	V5211	V5226	
Capture	V5203	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	
	V5204	-	-	+	-	+	-	-	-	-	-	-	-	-	-	-	-	
	V5205	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	V5206	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	V5207	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	V5212	+	-	-	-	-	+	+	-	-	-	-	-	-	-	-	-	
	V5213	+	+	+	+	+	-	-	-	-	-	-	-	+	+	+	-	
	V5215	+	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	
	V5216	+	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	
	V5217	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	V5218	-	-	-	-	-	+	-	-	-	-	-	+	-	-	-	-	
	V5220	-	-	-	-	-	-	-	-	-	-	+	-	-	+	-	-	
	V5231	-	-	-	-	-	-	-	+	+	-	-	-	-	-	-	-	
	V5233	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+++	
	V5234	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	
	V5235	-	-	-	-	-	+	+	-	-	+	-	-	-	-	-	-	
	V5236	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-	+++	
	W5201	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	
	V5225	-	+	+	-	-	+	-	-	-	+	-	-	-	-	-	-	
	V5227	+	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	
	V5222	-	-	-	+	-	+	-	-	-	-	+	-	-	-	-	-	
	V5223	-	-	+	-	-	+	-	-	-	-	+	-	-	-	-	-	