

BL0001

[Primary Antibody]

Blocker I Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Mouse

Isotype: IgG

Applications: ELISA (1:5000-10000)

Clonality: Monoclonal

CloneNo.: 7E2

Target: Blocker I

Purification: affinity purified by Protein A

Storage: PBS, pH7.6.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.