

**bsm-2025M**

**[ Primary Antibody ]**

## **RBC(G2H4) Mouse mAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Mouse

**Isotype:** IgG

**Clonality:** Monoclonal

**CloneNo.:** G2H4

**GeneID:** RBC

**Target:** RBC(G2H4)

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** Size : 50ul/100ul/200ul

0.01M PBS (pH7.4) with 0.02% Proclin300.

Size : 200ug (PBS only)

0.01M PBS

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Red blood cells are the most common type of blood cell and are the vertebrate body's principal means of delivering oxygen from the lungs or gills to body tissues via the blood.

**Applications:** **IHC-P** (1:100-500)

**IHC-F** (1:100-500)

**IF** (1:100-500)

**ELISA** (1:5000-10000)

**Reactivity:** (predicted: Human)