
VEGF165 Mouse mAb

Cat. No:	V2601
Name:	VEGF165 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	12H8K
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	>1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	Recombinant human VEGF165 Protein (bs-41613P)

PRODUCT SPECIFIC PUBLICATIONS

[IF=11.322] Yonghong Fan. et al. Construction of tissue-engineered vascular grafts with high patency by mimicking immune stealth and blocking TGF- β mediated endothelial-to-mesenchymal transition. COMPOS PART B-ENG. 2023 Feb;251:110487 IF ; Mouse . 10.1016/j.compositesb.2022.110487