
Ferritin Light Chain Rabbit pAb

Cat. No:	W2401
Name:	Ferritin Light Chain Rabbit pAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Polyclonal
Clone No:	
Host/Source:	Rabbit
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant human Ferritin Light Chain (bs-41189P)