

bs-6937R

[Primary Antibody]

PIG54 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 23244</p> <p>Target: PIG54</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SCC112/PIG54: 855-950/1337.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: The protein encoded by this gene binds to the cohesin complex and associates with chromatin through most of the cell cycle. The encoded protein may play a role in regulating sister chromatid cohesion during mitosis. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2010]</p>	<p>Isotype: IgG</p> <p>SWISS: Q29RF7</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:50-200) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Dog, Horse)</p> <p>Predicted MW.: 147 kDa</p> <p>Subcellular Location: Nucleus</p>
---	--	--