## bs-6937R

## [ Primary Antibody ]

## BIOSS ANTIBODIES

## PIG54 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**GenelD:** 23244 **SWISS:** Q29RF7

Target: PIG54

Immunogen: KLH conjugated synthetic peptide derived from human

SCC112/PIG54: 855-950/1337.

**Purification:** affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

**Background:** The protein encoded by this gene binds to the cohesin complex

and associates with chromatin through most of the cell cycle. The encoded protein may play a role in regulating sister chromatid cohesion during mitosis. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May

2010]

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Rabbit, Pig, Sheep,

Dog, Horse)

Predicted MW.: 147 kDa

Subcellular Nucleus