

bs-19167R

[Primary Antibody]

MYO1B Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal GeneID: 4430 Target: MYO1B Immunogen: KLH conjugated synthetic peptide derived from human MYO1B: 781-880/1136. Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Background: Myosin Ib is a motor protein that is involved in cell migration, neurite outgrowth and vesicular transport. In multivesicular endosomes, Myosin Ib has been implicated in protein cargo traffic control.	Isotype: IgG SWISS: O43795	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse, Chinese Hamster) Predicted MW.: 132 kDa Subcellular Location: Cytoplasm
--	--	---