

bs-10654R

[Primary Antibody]

Duck Hepatitis Virus VP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Duck Hepatitis Virus VP1

Immunogen: KLH conjugated synthetic peptide derived from Duck Hepatitis Virus VP1: 101-200/240.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Duck hepatitis is caused by the enteroviruses DHV-1 and DHV-3. It is a fatal disease of ducklings causing opisthotonus and hepatitis. Duck hepatitis virus, DHV, type 1 (DHV-1), type 1a (DHV-1a), type 2 (DHV-2), and type 3 (DHV-3) are the four types of viruses that can cause disease in ducks. The pathology and pathogenicity between duck hepatitis virus type 1 and type 1a is, in essence, the same as one another.

Applications: **WB** (1:500-2000)

IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: DHAV)

Predicted MW.: 26 kDa