

---

## DNER Rabbit pAb

Catalog Number: bs-12115R

Target Protein: DNER

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Sheep, Cow, Dog)

Predicted MW: 75 kDa

Entrez Gene: 92737

Swiss Prot: Q8NFT8

Source: KLH conjugated synthetic peptide derived from human DNER: 185-300/737.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Activator of the NOTCH1 pathway. May mediate neuron-glia interaction during astrocytogenesis.