
phospho-SHB (Tyr246) Rabbit pAb

Catalog Number: bs-19758R

Target Protein: phospho-SHB (Tyr246)

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Cow, Horse)

Predicted MW: 55 kDa

Entrez Gene: 6461

Swiss Prot: Q15464

Source: KLH conjugated synthesised phosphopeptide derived from human SHB around the phosphorylation site of Tyr246: Y(p-S)DP.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.