

**bs-18665R****[ Primary Antibody ]****MAP3K6/ASK2 Rabbit pAb****BioSS**  
**ANTIBODIES**

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**— DATASHEET —**

|                                                                                                                                                                                                                                                                                                                                                                |                      |                                                                                                                                                                                                                                                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host:</b> Rabbit                                                                                                                                                                                                                                                                                                                                            | <b>Isotype:</b> IgG  | <b>Applications:</b> IHC-P (1:100-500)<br>IHC-F (1:100-500)<br>IF (1:100-500)<br><br><b>Reactivity:</b> Mouse (predicted: Human, Rat, Dog, Horse)<br><br><b>Predicted MW.:</b> 143 kDa<br><br><b>Subcellular Location:</b> Cell membrane ,Cytoplasm |
| <b>Clonality:</b> Polyclonal                                                                                                                                                                                                                                                                                                                                   |                      |                                                                                                                                                                                                                                                     |
| <b>GeneID:</b> 9064                                                                                                                                                                                                                                                                                                                                            | <b>SWISS:</b> O95382 |                                                                                                                                                                                                                                                     |
| <b>Target:</b> MAP3K6/ASK2                                                                                                                                                                                                                                                                                                                                     |                      |                                                                                                                                                                                                                                                     |
| <b>Immunogen:</b> KLH conjugated synthetic peptide derived from human MAP3K6/ASK2: 101-200/1288.                                                                                                                                                                                                                                                               |                      |                                                                                                                                                                                                                                                     |
| <b>Purification:</b> affinity purified by Protein A                                                                                                                                                                                                                                                                                                            |                      |                                                                                                                                                                                                                                                     |
| <b>Concentration:</b> 1mg/ml                                                                                                                                                                                                                                                                                                                                   |                      |                                                                                                                                                                                                                                                     |
| <b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                                                                                                                                                                           |                      |                                                                                                                                                                                                                                                     |
| <b>Background:</b> This gene encodes a serine/threonine protein kinase that forms a component of protein kinase-mediated signal transduction cascades. The encoded kinase participates in the regulation of vascular endothelial growth factor (VEGF) expression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014] |                      |                                                                                                                                                                                                                                                     |

**— VALIDATION IMAGES —**

Paraformaldehyde-fixed, paraffin embedded (Mouse brain); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (MAP3K6) Polyclonal Antibody, Unconjugated (bs-18665R) at 1:400 overnight at 4°C, followed by operating according to SP Kit(Rabbit) (sp-0023) instructions and DAB staining.