

bs-14351R

[Primary Antibody]

DMRTD2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Dog) Predicted MW.: 39 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 63946	SWISS: Q8IXT2	
Target: DMRTD2		
Immunogen: KLH conjugated synthetic peptide derived from human DMRTD2: 261-367/367.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=4.183]** Li TT et al. Cloning, Molecular Characterization and Expression Patterns of DMRTC2 Implicated in Germ Cell Development of Male Tibetan Sheep. Int J Mol Sci. 2020 Apr 1;21(7):2448. WB,IF ;Sheep. 32244802