

bs-14609R

[Primary Antibody]

EPG5 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 57724</p> <p>Target: EPG5</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human EPG5: 5-100/2579.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X. [provided by RefSeq, Mar 2009]</p>	<p>Isotype: IgG</p> <p>SWISS: Q9HCE0</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Dog, Horse)</p> <p>Predicted MW.: 292 kDa</p> <p>Subcellular Location: Cell membrane ,Nucleus</p>
--	--	--