

---

## NSUN3 Rabbit pAb

Catalog Number: bs-19478R

Target Protein: NSUN3

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse)

Predicted MW: 38 kDa

Entrez Gene: 63899

Swiss Prot: Q9H649

Source: KLH conjugated synthetic peptide derived from human NSUN3: 1-100/340.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.