

bs-19478R

[Primary Antibody]

NSUN3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|---|--|---|
| <p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 63899</p> <p>Target: NSUN3</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human NSUN3: 1-100/340.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> | <p>Isotype: IgG</p> <p>SWISS: Q9H649</p> | <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse)</p> <p>Predicted MW.: 38 kDa</p> <p>Subcellular Location: Nucleus</p> |
|---|--|---|