

---

## Transketolase (CT) Rabbit pAb

Catalog Number: bs-4514R

Target Protein: Transketolase (CT)

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:1000-5000)

Reactivity: (predicted:Apple)

Predicted MW: 82 kDa

Source: KLH conjugated synthetic peptide derived from Apple Transketolase: 621-720/747.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.