bs-9703R

[Primary Antibody]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

S100A18 Rabbit pAb

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 388697 SWISS: Q86YZ3

Target: S100A18

Immunogen: KLH conjugated synthetic peptide derived from human

Hornerin/S100A18: 1-100/2850.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: Hornerin has similar structural features to those of profilaggrin, an

essential protein for keratinization of epidermal tissues. Hornerin

may play a role in cornification.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) **IF** (1:100-500) **ELISA** (1:5000-10000)

Reactivity: (predicted: Human)

Predicted MW.: ^{282 kDa}

Subcellular Cytoplasm