

---

## S100A18 Rabbit pAb

Catalog Number: bs-9703R

Target Protein: S100A18

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: **IHC-P** (1:100-500), **IHC-F** (1:100-500), **IF** (1:50-200), **ELISA** (1:5000-10000)

Reactivity: (predicted:&nbsp;Human)

Predicted MW: 282 kDa

Subcellular Cytoplasm

Locations:

Entrez Gene: 388697

Swiss Prot: Q86YZ3

Source: KLH conjugated synthetic peptide derived from human Hornerin/S100A18: 1-100/2850.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Hornerin has similar structural features to those of profilaggrin, an essential protein for keratinization of epidermal tissues. Hornerin may play a role in cornification.