

bs-9703R

[Primary Antibody]

S100A18 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|--|--|
| <p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 388697</p> <p>Target: S100A18</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human Hornerin/S100A18: 1-100/2850.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Hornerin has similar structural features to those of profilaggrin, an essential protein for keratinization of epidermal tissues. Hornerin may play a role in cornification.</p> | <p>Isotype: IgG</p> <p>SWISS: Q86YZ3</p> | <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:50-200) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 282 kDa</p> <p>Subcellular Location: Cytoplasm</p> |
|--|--|--|