

bs-2328R

[Primary Antibody]

HBV pre S1 protein Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>Target: HBV pre S1 protein</p> <p>Immunogen: KLH conjugated synthetic peptide derived from Hepatitis B Virus Surface Antigen: 2-80/389.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background:</p>	<p>Isotype: IgG</p> <p>Applications: ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Hepatitis B virus(HBVD2))</p> <p>Predicted MW.: 31/44 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
---	--