

---

## HBV pre S1 protein Rabbit pAb

Catalog Number: bs-2328R

Target Protein: HBV pre S1 protein

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:Hepatitis B virus(HBVD2))

Predicted MW: 31/44 kDa

Source: KLH conjugated synthetic peptide derived from Hepatitis B Virus Surface Antigen: 2-80/389.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: