

bs-9942R**[Primary Antibody]****C11orf73 Rabbit pAb****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000)
Clonality: Polyclonal		Reactivity: Human, Rat (predicted: Mouse, Rabbit, Pig, Sheep, Cow, Chicken, Dog)
GeneID: 51501	SWISS: Q53FT3	Predicted MW.: 22 kDa
Target: C11orf73		Subcellular Location: Cytoplasm ,Nucleus
Immunogen: KLH conjugated synthetic peptide derived from human C11orf73: 21-120/197.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Required for organization and/or function of the secretory apparatus in Clara cells in lung.		

— VALIDATION IMAGES —

25 ug total protein per lane of various lysates (see on figure) probed with C11orf73 polyclonal antibody, unconjugated (bs-9942R) at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.